

Peroxichrom™, Peroxichrom Excel™

TMB peroxidase substrate solution (1-component) for ELISA

Cat. # :

D5015-100 (Peroxichrom™ : 100ml)


D5015-400 (Peroxichrom™ : 400ml)

D5016-100 (peroxichrom Excel™ : 100ml)

D5016-400 (peroxichrom Excel™ : 400ml)

Highlights of TMB substrate kit

- ▶ Ready to use single component
- ▶ Highest sensitivity
- ▶ Sufficient dynamic range
- ▶ Easy to use
- ▶ Noncarcinogenic
- ▶ No DMF or DMSO present in the reagent
- ▶ Stable at RT for 1 month (Stable at 4°C for 1 year)
- ▶ Easy to transport



Description:

TMB substrate(3,3'.5,5'-tetramethylbenzidine) is a chromogen that yields a deep blue color (maximum absorbance at 605nm or 650nm) when oxidized with hydrogen peroxide(catalyzed by HRP). The color then changes to yellow with the addition of 2N H₂SO₄ with maximum absorbance at 450 nm. Our TMB Substrates(Peroxichrom™, Peroxichrom Excel™) are one-component substrates that require no preparation before using, stable and sensitivity. Also Peroxichrom™ and Peroxichrom Excel™ contain no organic solvents such as DMF, DMSO, and methanol so there is no issue of safety with user.

Precautions:

TMB substrate is sensitive to contamination from a variety of oxidizing agents. Avoid prolonged exposure to light, contact with metal or air. We recommend using TMB substrate by pouring out required amount into a reservoir and do not return excess TMB to provided bottle.

Storage:

Stable at RT for 1 month (Stable at 4°C for 1 year)

Procedure:

1. Warm to room temperature prior to use if you store TMB solution at 4°C
2. Dispense 1-Component appropriate TMB solution or 100ul into each well
3. After sufficient color development (5-10minutes at room temperature or at 37°C)
Add 50ul Stop Solution (2N H₂SO₄ or 1M H₃PO₄) to each well.
4. Read plates at 450nm